GYROMAX™ 727 ORBITAL INCUBATOR SHAKER

Incubate and Shake it to the Max...with GYROMAX!

The GYROMAXTM 727 Orbital Incubator Shaker is compact and occupies only 20"W x



27.5"D of space. With a height of 30.5", it can fit on or under a workbench. The 727 has an adjustable shelf (one is included) for simultaneous incubation with or without shaking. Up to 6 shelves can fit in the chamber. The chamber holds a 15"W x 15"D platform and accommodates flasks up to 4 liters. The double-wall construction, together with 1.25" of insulation and a triple-pane glass window, minimizes heat loss. This construction, together with an efficient circulation system, ensures unrivaled optimum temperature uniformity throughout the chamber. The temperature is microprocessor controlled monitored accurately by a Pt-100 RTD. It is set digitally and reproducibly from a membrane keypad and displayed by an LED to 0.1°C. The working temperature is from about 7°C above ambient to 80°C. An independent over-temperature thermostat protects valuable samples overheating. The from timer programmable up to 99 hours and 59 minutes or for continuous operation.

The GYROMAX Incubator Shaker features a triple-eccentric drive with counter-balancing weights for smooth and quiet shaking of heavy loads without the need for adjustments. A large, quiet fan ensures air temperature uniformity throughout the stainless steel chamber. A brushless, maintenance-free, drive motor generates digitally displayed speeds from 40 to 400 RPM. The standard shaking orbit of 25 mm provides a higher oxygen transfer rate (OTR) compared to shakers with a 19-mm orbit because OTR is directly proportional to orbit size. The shaking stops when the door is opened. The chamber lamp is turned on and off with a switch on the control panel. The shaker can accommodate Erlenmeyer flasks up to four liters.

Optional Accessories Include: platforms, clamps, spring racks, tube racks, and plate holders. Platform holes are drilled to accommodate the different clamps purchased. The spring rack is designed for holding tubes in a horizontal, angled, or an upright position. It also can be used to hold flasks up to 1-liter without clamps. The rack consists of a 2-tier grid of springs. The distance between spring racks is adjustable to accommodate tubes and flasks of different heights. They are mounted onto the platform by four screws.

GYROMAX™ 727 BENCHTOP INCUBATOR SHAKER

SPECIFICATIONS	727					
Temperature Range	Ambient +7°C to 80°C					
Temperature Control	PID microprocessor controller, Pt-100 RTD sensor					
Accuracy/Uniformity	±0.1°C/±0.5°C at 37.0°C					
Temperature Display						
Speed Range	40 to 400 RPM with standard 1" (25mm) orbit					
Speed Display	Digital					
Orbit	25mm (1") orbit provides higher oxygen transfer rate than 19mm					
Timer	Programmable up to 99 hours 59 minutes or for continuous operation					
Insulation	1.25" double-wall housing w/ polyurethane insulation minimizes heat					
	loss					
Safety Devices	An independent thermostat prevents overheating. Shaking stops					
	when door is opened. Auto-restart after power failure.					
Stacking	Two can be stacked to save space.					
Drive	Dynamically balanced triple-eccentric drive with counter-balancing					
	weights provides smooth and quiet shaking without the need to make					
Barton	inconvenient adjustments as weight of the load changes					
Motor	Maintenance-free brushless motor					
Shelf	One included, five additional shelves optional					
Viewing Window	Triple-pane 1.25" (30mm) thick					
Spring Rack	Accommodates an accessory spring rack that can shake flasks up to					
	1-L without the use of clips or clamps as well as 15-mL or 50-mL					
	tubes placed vertically or at 45 ^o to the platform					
Chamber Material	Stainless steel					
Platform Dimensions	15"W x 15"D (380mm x 380mm)					
Chamber Dimensions	17.7"W x 17.9"D x 17.3"H (450mm x 455mm x 440mm)					
Overall Dimensions	20"W x 27.5"D x 30.5"H (510mm x 700mm x 770mm)					
Power Source	120V (or 220V), 50/60Hz, 6A					
Warranty	2 years					

Flask Capacity

Flask Size (mL)	25	50	125	250	500	1000	2000	2800	4000
Number of Flasks	49	36	25	16	9	5	4	2	1

Accessories:

Holders for Tube Racks & Plates

Full Size Spring Rack









Flask Clamps





Amerex Instruments, Inc. 3951C Industrial Way Concord, California 94520 925.299.0743 / Fax: 925.299.0745 www.amerexinst.com marketing@amerexinst.com